

Product datasheet for RC219886

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

CD30 (TNFRSF8) (NM_152942) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: CD30 (TNFRSF8) (NM_152942) Human Tagged ORF Clone

Tag: Myc-DDK
Symbol: TNFRSF8

Synonyms: CD30; D1S166E; Ki-1

Vector:pCMV6-Entry (PS100001)E. coli Selection:Kanamycin (25 ug/mL)

Cell Selection: Neomycin

ORF Nucleotide >RC219886 representing NM_152942

Sequence: Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC

GCCGCGATCGCC

ATGAGCCAGCCACTGATGGAGACCTGCCACAGCGTGGGGGCAGCCTACCTGGAGAGCCTGCCGCTGCAGG ATGCCAGCCCGGCCGGGGGCCCCTCGTCCCCCAGGGACCTTCCTGAGCCCCGGGTGTCCACCGGAGCACAC CAATAACAAGATTGAGAAAATCTACATCATGAAGGCTGACACCGTGATCGTGGGGACCGTGAAGGCTGAG CTGCCGGAGGGCCGGGGCCTGGCGGGGCCAGCAGAGCCCGAGTTGGAGGAGGAGCTGGAGGCGGACCATA CCCCCCACTACCCCGAGCAGGAGACAGAACCGCCTCTGGGCAGCTGCAGCGATGTCATGCTCTCAGTGGA

AGAGGAAGGAAAGAACCCCTTGCCCACAGCTGCCTCTGGAAAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT

ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC219886 representing NM_152942

Red=Cloning site Green=Tags(s)

MSQPLMETCHSVGAAYLESLPLQDASPAGGPSSPRDLPEPRVSTEHTNNKIEKIYIMKADTVIVGTVKAE

LPEGRGLAGPAEPELEEELEADHTPHYPEQETEPPLGSCSDVMLSVEEEGKEDPLPTAASGK

TRTRPLEQKLISEEDLAANDILDYKDDDDK**V**


Restriction Sites: Sgfl-Mlul






Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

Plasmid Map:



ACCN: NM_152942

ORF Size: 396 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info



CD30 (TNFRSF8) (NM_152942) Human Tagged ORF Clone - RC219886

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 152942.2, NP 694421.1</u>

RefSeq Size: 2361 bp
RefSeq ORF: 398 bp
Locus ID: 943
Cytogenetics: 1p36.22

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Stem cell - Pluripotency, Transmembrane

Protein Pathways: Cytokine-cytokine receptor interaction

MW: 13.9 kDa


Gene Summary: The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor

is expressed by activated, but not by resting, T and B cells. TRAF2 and TRAF5 can interact with this receptor, and mediate the signal transduction that leads to the activation of NF-kappaB. This receptor is a positive regulator of apoptosis, and also has been shown to limit the proliferative potential of autoreactive CD8 effector T cells and protect the body against autoimmunity. Two alternatively spliced transcript variants of this gene encoding distinct

isoforms have been reported. [provided by RefSeq, Jul 2008]



Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY TNFRSF8 (Cat# RC219886, Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TNFRSF8(Cat# [TA590389]).