

Product datasheet for **RC219888A1V**

Human ERK1 (MAPK3) (NM_001040056) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ERK1 (MAPK3) (NM_001040056) AAV Particle
Tag:	Myc-DDK
Symbol:	ERK1
Synonyms:	ERK-1; ERK1; ERT2; HS44KDAP; HUMKER1A; p44-ERK1; p44-MAPK; P44ERK1; P44MAPK; PRKM3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC219888 representing NM_001040056 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGGGCGGGCGGCTCAGGGGGCGGGGGCGGGGAGCCCCGTAGAACCAGGGGGTTCGGCCCGGGG
TCCCGGGGAGGTGGAGATGGTGAAGGGGCGCCGTTTCGACGTGGGCCCGCCTACACGCAGTTGCAGTA
CATCGGCGAGGGCGCGTACGGCATGGTCACTCGGCTATGACCACGTGCGCAAGACTCGCGTGGCCATC
AAGAAGATCAGCCCTTCGAACATCAGACCTACTGCCAGCGCACGCTCCGGGAGATCCAGATCCTGCTGC
GCTTCGCCCATGAGAATGTCATCGGCATCCGAGACATTCGCGGGCGTCCACCCTGGAAGCCATGAGAGA
TGTCTACATTGTGCAGGACCTGATGGAGACTGACCTGTACAAGTTGCTGAAAAGCCAGCAGCTGAGCAAT
GACCATATCTGCTACTTCTCTACCAGATCCTGCGGGGCTCAAGTACATCCACTCCGCCAACGTGCTCC
ACCGAGATCTAAAGCCCTCCAACCTGCTCATCAACACCACCTGCGACCTTAAGATTTGTGATTCGGCCT
GGCCCGGATTGCCGATCCTGAGCATGACCACACCGGCTTCTGACGGAGTATGTGGCTACGCGTGGTAC
CGGGCCCCAGAGATCATGCTGAACCTCAAGGGCTATACCAAGTCCATCGACATCTGGTCTGTGGGCTGCA
TTCTGGCTGAGATGCTCTTAACCGGCCATCTCCCTGGCAAGCACTACCTGGATCAGCTCAACCACAT
TCTGGCATCCTGGGCTCCCCATCCCAGGAGGACCTGAATTGTATCATCAACATGAAGGCCCGAACTAC
CTACAGTCTCTGCCCTCAAGACCAAGTGGCTTGGGCCAAGCTTTTCCCAAGTCAAGTCCAAAGCCC
TTGACCTGCTGGACCGGATGTTAACCTTTAACCCAATAAACGGATCACAGTGGAGGAAGCGCTGGCTCA
CCCCTACCTGGAGCAGTACTATGACCCGACGGATGAGGTGGGCCAGTCCCAGCAGCAGTGGGGCTGGG
GCAGGGGAGCAGGGGGGCACG

ACGCGTACGCGGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC219888 representing NM_001040056 Red=Cloning site Green=Tags(s) MAAAAAQGGGGGPRRTEGVGPGVPEVEMVKGQPFVDVGPRYTQLQYIGEGAYGMVSSAYDHVRKTRVAI KKISPFHQTYCQRTLREIQILLRFRHENVIGIRDILRASTLEAMRDVYIVQDLMETDL YKLLKSQQLSN DHICYFLYQILRGLKYIHSANVLRDLKPSNLLINTTCDLKICDFGLARIADPEHDHTGFLTEYVATRWY RAPEIMLNSKGYTKSIDIWSVGCILAEMLSNRPIFPKGHYLDQLNHILGILGSPSQEDLNCIINMKARNY LQSLPSKTKVAVAKLFPKSDSKALDLLDRMLTFNPNKRITVEEALAHPLYEQYYDPTDEVGQSPAAVGLG AGEQGGT TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001040056
ORF Size:	1071 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001040056.1
RefSeq Size:	2005 bp
RefSeq ORF:	1074 bp
Locus ID:	5595
UniProt ID:	P27361
Cytogenetics:	16p11.2
MW:	39.9 kDa