

Product datasheet for **RC219889A1V**

Human CXorf27 (HYPM) (NM_012274) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CXorf27 (HYPM) (NM_012274) AAV Particle
Tag: Myc-DDK
Symbol: CXorf27
Synonyms: CXorf27; HIP17; HYPM
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC219889 representing NM_012274
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGTCAGAGAAAAAGAACTGTAAGAACTCGTCTACAAATAACAACCAGACTCAAGACCCTTCTAGAAATG
 AGCTACAGGTCCTTAGGAGCTTCGTGGACCGCGTTGTGCAAGATGAACGAGACGTCCAAAGCCAGAGTTC
 CTCCACAATAAATACCCTCCTTACCTTGCTCGATTGCCTTGCTGACTACATCATGGAGCGGGTAGGCTTG
 GAGGCCAGCAACAATGGCAGTATGCGCAACACTTCACAAGATAGGGAGAGAGAAGTGGACAACAACCGTG
 AGCCCCACAGCGCTGAGAGTGATGTGACTCGCTTTTGTGTTGATGAGATGCCCAAATCCAGGAAGAATGA
 C

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC219889 representing NM_012274
 Red=Cloning site Green=Tags(s)

MSEKKKNKNSSTNNNQTDPSRNLQVPRSFVDRVVQDERDVQSQSSSTINTLLTLLDCLADYIMERVGL
 EASNNGSMRNTSQDREREVDNNREPHSAESDVTRFLFDEMPKSRKND

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_012274



[View online »](#)

ORF Size:	351 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_012274.1</u>
RefSeq Size:	517 bp
RefSeq ORF:	354 bp
Locus ID:	25763
UniProt ID:	<u>O75409</u>
Cytogenetics:	Xp11.4
MW:	13.3 kDa