

## Product datasheet for **RC219893L4V**

### SCML1 (NM\_001037536) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SCML1 (NM_001037536) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SCML1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001037536
ORF Size:	906 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219893).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001037536.2</a> , <a href="#">NP_001032625.1</a>
RefSeq Size:	2701 bp
RefSeq ORF:	627 bp
Locus ID:	6322
UniProt ID:	<a href="#">Q9UN30</a>
Cytogenetics:	Xp22.13
Protein Families:	Transcription Factors
MW:	34.2 kDa



[View online »](#)

**Gene Summary:**

Putative Polycomb group (PcG) protein. PcG proteins act by forming multiprotein complexes, which are required to maintain the transcriptionally repressive state of homeotic genes throughout development. May be involved in spermatogenesis during sexual maturation (By similarity).[UniProtKB/Swiss-Prot Function]