

## Product datasheet for **RC219897L3V**

### ZFYVE27 (NM\_001002261) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ZFYVE27 (NM_001002261) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZFYVE27
Synonyms:	PROTRUDIN; SPG33
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001002261
ORF Size:	1248 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219897).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001002261.1</a> , <a href="#">NP_001002261.1</a>
RefSeq Size:	3045 bp
RefSeq ORF:	1251 bp
Locus ID:	118813
UniProt ID:	<a href="#">Q5T4F4</a>
Cytogenetics:	10q24.2
Protein Families:	Transmembrane
MW:	46.2 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein with several transmembrane domains, a Rab11-binding domain and a lipid-binding FYVE finger domain. The encoded protein appears to promote neurite formation. A mutation in this gene has been reported to be associated with hereditary spastic paraplegia, however the pathogenicity of the mutation, which may simply represent a polymorphism, is unclear. [provided by RefSeq, Mar 2010]