

Product datasheet for **RC219907L4V**

C6orf89 (NM_152734) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C6orf89 (NM_152734) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C6orf89
Synonyms:	BRAP; PS1TP5TP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_152734
ORF Size:	1062 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219907).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_152734.2
RefSeq Size:	2271 bp
RefSeq ORF:	1065 bp
Locus ID:	221477
UniProt ID:	Q6UWU4
Cytogenetics:	6p21.2
Protein Families:	Transmembrane
MW:	40.6 kDa



[View online »](#)

Gene Summary:

Exhibits histone deacetylase (HDAC) enhancer properties (PubMed:23460338). May play a role in cell cycle progression and wound repair of bronchial epithelial cells (PubMed:21857995). [UniProtKB/Swiss-Prot Function]