

Product datasheet for RC219915L4V

SAMD11 (NM_152486) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SAMD11 (NM_152486) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SAMD11
Synonyms:	MRS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_152486
ORF Size:	2043 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219915).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_152486.2
RefSeq Size:	2575 bp
RefSeq ORF:	2046 bp
Locus ID:	148398
UniProt ID:	Q96NU1
Cytogenetics:	1p36.33
MW:	72.5 kDa
Gene Summary:	May play a role in photoreceptor development.[UniProtKB/Swiss-Prot Function]


[View online »](#)