

Product datasheet for **RC219918L4V**

p90 Autoantigen (KIAA1524) (NM_020890) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	p90 Autoantigen (KIAA1524) (NM_020890) Human Tagged ORF Clone Lentiviral Particle
Symbol:	p90 Autoantigen
Synonyms:	KIAA1524; p90
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020890
ORF Size:	2715 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219918).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020890.1
RefSeq Size:	3395 bp
RefSeq ORF:	2718 bp
Locus ID:	57650
UniProt ID:	Q8TCG1
Cytogenetics:	3q13.13
Protein Families:	Druggable Genome
MW:	102 kDa


[View online »](#)

Gene Summary:

Oncoprotein that inhibits PP2A and stabilizes MYC in human malignancies. Promotes anchorage-independent cell growth and tumor formation.[UniProtKB/Swiss-Prot Function]