

Product datasheet for RC219955L4

USP35 (NM_020798) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	USP35 (NM_020798) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	USP35
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219955).
Restriction Sites:	Sgfl-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_020798
ORF Size:	3054 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_020798.1
RefSeq Size:	4197 bp
RefSeq ORF:	3057 bp
Locus ID:	57558
UniProt ID:	Q9P2H5
Cytogenetics:	11q14.1
Protein Families:	Druggable Genome, Protease
MW:	113.2 kDa
Gene Summary:	This gene encodes a member of the peptidase C19 family of ubiquitin-specific proteases. These deubiquitinating enzymes (DUBs) catalyze the removal of ubiquitin proteins from other proteins. The encoded protein associates with polarized mitochondria and has been shown to inhibit NF-kappa B activation and delay PARK2-mediated degradation of mitochondria. Expression of this gene is upregulated by the let-7a microRNA and reduced expression has been observed in human tumor tissues. [provided by RefSeq, Jul 2017]

Product images:



Circular map for RC219955L4