

Product datasheet for **RC220043L3V**

HS6ST2 (NM_001077188) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HS6ST2 (NM_001077188) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HS6ST2
Synonyms:	MRXSPM
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001077188
ORF Size:	1935 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220043).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001077188.1
RefSeq Size:	4576 bp
RefSeq ORF:	1938 bp
Locus ID:	90161
UniProt ID:	Q96MM7
Cytogenetics:	Xq26.2
Protein Families:	Transmembrane
Protein Pathways:	Heparan sulfate biosynthesis



[View online »](#)

MW: 73.2 kDa

Gene Summary: Heparan sulfate proteoglycans are ubiquitous components of the cell surface, extracellular matrix, and basement membranes, and interact with various ligands to influence cell growth, differentiation, adhesion, and migration. This gene encodes a member of the heparan sulfate (HS) sulfotransferase gene family, which catalyze the transfer of sulfate to HS. Different family members and isoforms are thought to synthesize heparan sulfates with tissue-specific structures and functions. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]