

Product datasheet for RC220120L3

KIF24 (NM_194313) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KIF24 (NM_194313) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	KIF24
Synonyms:	bA571F15.4; C9orf48
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220120).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_194313
ORF Size:	4104 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_194313.2
RefSeq Size:	6237 bp
RefSeq ORF:	4107 bp
Locus ID:	347240
UniProt ID:	Q5T7B8
Cytogenetics:	9p13.3
Protein Families:	Druggable Genome
MW:	151.7 kDa
Gene Summary:	This gene encodes a member of the kinesin superfamily of microtubule-based motor proteins which are involved in the intracellular transport of membranous organelles, protein complexes, and mRNAs. They also play critical roles in mitosis, morphogenesis, and signal transduction. The encoded protein contains an N-terminal sterile alpha motif (SAM) domain and an ATP-binding kinesin motor domain. It binds centriolar coiled coil protein 110 and centrosomal protein 97 and localizes to the mother centriole to regulate ciliogenesis by controlling microtubule polymerization. [provided by RefSeq, Mar 2017]

Product images:



Circular map for RC220120L3