

Product datasheet for **RC220127A1V**

Human PLSCR5 (NM_001085420) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human PLSCR5 (NM_001085420) AAV Particle
Tag: Myc-DDK
Symbol: PLSCR5
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC220127 representing NM_001085420
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCTCTAAAGATGCCAGAACCAAGAAGAGGTCTGCCTGGTTTTCTTCTGGAGCTCCAGACCCAG
ACCAAGCCTTCTGCCTCTTCCAATCCAGGGAACCAAGCATGGCAGCTGAGTCTCCCTCTGCCAAGCAG
TTTCTGCCAACAGTCAGTCTCCCTCCTGGTCTAGAATATTTAAGCCAGTTAGACCTGATAATTATACAC
CAGCAGGTGGAGCTGCTTGAATGATACTTGGTACTGAGACCTCCAACAAATATGAGATTAACAAACAGCT
TGGGACAAAGAATTTACTTTGCAGTGGAGGAAAGCATCTGCTTCAATCGTACTTTCTGTTCCACTCTGCG
ATCTTGACCCCTGAGGATCACAGATAACTCAGGTGAGAGGTCATTACAGTGAACAGGCCCTTGAGATGT
AACAGCTGCTGGTGCCCTTGCTACCTACAAGAGTTAGAAATCCAAGCCCCCTCTGGTACTATAGTTGGTT
ACGTTACGCAGAAGTGGGACCCCTTTCTGCCTAAATTCACAATCCAAATGCAAACAAAGAAGATATTTT
GAAAAATTGTTGGTCCTTGTGTGACATGTGGCTGTTTTGGCGATGTGGATTTTGAGGTGAAACCATTAAT
GAAAAGCTTACAATTGGAAGATTTCAAAGTACTGGTCAGGATTTGTAATGATGTCTTCACAAATGCTG
ACAATTTCCGAATTCATGTTCTGTCAGATCTAGATGTAAACAGTCAAAGCAGCAATGATCGGTGCCTGTTT
TCTCTTTGATTTTATGTTCTTTGAACATTCAGTGGCTGGATTA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC220127 representing NM_001085420 Red=Cloning site Green=Tags(s) MASKDAQNQRRGLPGFLPGAPDPDQSLPASSNPGNQAWQLSLPLPSSFLPTVSLPPGLELYLSQLDLIIH QQVELLGMILGTETSNKYEIKNSLGQRIYFAVEESICFNRTFCSTLRCTLRITDMSGREVITVNRPLRC NSWCPCYLQELEIQAPPGTIVGYVTQKWDPFLPKFTIQNANKEDILKIVGPCVTCGCFGDVDFEVKTIN EKLTI GKISKYWSGFVNDVFTNADNFGIHVPADLDVTVKAAMIGACFLFDFMFFEHSAGL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001085420
ORF Size:	813 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001085420.1, NP_001078889.1</u>
RefSeq Size:	1969 bp
RefSeq ORF:	816 bp
Locus ID:	389158
UniProt ID:	<u>A0PG75</u>
Cytogenetics:	3q24
MW:	29.8 kDa