

Product datasheet for **RC220128L3V**

SLK (NM_014720) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLK (NM_014720) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLK
Synonyms:	bA16H23.1; LOSK; se20-9; STK2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_014720
ORF Size:	3705 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220128).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014720.2
RefSeq Size:	5629 bp
RefSeq ORF:	3708 bp
Locus ID:	9748
UniProt ID:	Q9H2G2
Cytogenetics:	10q24.33-q25.1
Domains:	pkinase
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

Protein Pathways: Oocyte meiosis

MW: 142.5 kDa

Gene Summary: Mediates apoptosis and actin stress fiber dissolution.[UniProtKB/Swiss-Prot Function]