

## Product datasheet for **RC220131L3V**

### Folate Binding Protein (FOLR1) (NM\_016725) Human Tagged ORF Clone Lentiviral Particle

#### Product data:


Product Type:	Lentiviral Particles
Product Name:	Folate Binding Protein (FOLR1) (NM_016725) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Folate Binding Protein
Synonyms:	FBP; FOLR; FRalpha; NCFTD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016725
ORF Size:	771 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220131).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_016725.1</a>
RefSeq Size:	1104 bp
RefSeq ORF:	774 bp
Locus ID:	2348
UniProt ID:	<a href="#">P15328</a>
Cytogenetics:	11q13.4
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
MW:	29.82 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a member of the folate receptor family. Members of this gene family bind folic acid and its reduced derivatives, and transport 5-methyltetrahydrofolate into cells. This gene product is a secreted protein that either anchors to membranes via a glycosyl-phosphatidylinositol linkage or exists in a soluble form. Mutations in this gene have been associated with neurodegeneration due to cerebral folate transport deficiency. Due to the presence of two promoters, multiple transcription start sites, and alternative splicing, multiple transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Oct 2009]

**Product images:**

[RC220131L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC220131L3V particle to overexpress human FOLR1-Myc-DDK fusion protein.