

Product datasheet for **RC220159L3V**

ADCK2 (NM_052853) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ADCK2 (NM_052853) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ADCK2
Synonyms:	AARF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_052853
ORF Size:	1878 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220159).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_052853.3
RefSeq Size:	2376 bp
RefSeq ORF:	1881 bp
Locus ID:	90956
UniProt ID:	Q7Z695
Cytogenetics:	7q34
Domains:	ABC1
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

MW: 68.8 kDa

Gene Summary: The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr).[UniProtKB/Swiss-Prot Function]