

## Product datasheet for **RC220213**

### NRG1 (NM\_013957) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	NRG1 (NM_013957) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NRG1
Synonyms:	ARIA; GGF; GGF2; HGL; HRG; HRG1; HRGA; MST131; MSTP131; NDF; NRG1-IT2; SMDF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RC220213 representing NM\_013957  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGTCGAGCGCAAAGAAGGCAGAGGCAAAGGGAAGGCAAGAAGAAGGAGCGAGGCTCCGCAAGAAGC  
 CGGAGTCCGCGGGCAGCCAGAGCCAGCCTTGCCCTCCCGATTGAAAGAGATGAAAAGCCAGGAATC  
 GGCTGCAGGTTCCAACTAGTCCTTCGGTGTGAAACCAGTTCTGAATACTCCTCTCAGATTCAAGTGG  
 TTCAAGAATGGGAATGAATTGAATCGAAAAACAAACCACAAAATATCAAGATACAAAAAGCCAGGGA  
 AGTCAGAACTTCGATTAACAAAGCATCACTGGCTGATTCTGGAGAGTATATGTGCAAAGTGATCAGCAA  
 ATTAGGAAATGACAGTGCCTCTGCCAATATCACCATCGTGAATCAAACGAGATCATCACTGGTATGCCA  
 GCCTCACTGAAGGAGCATATGTGTCTCAGAGTCTCCATTAGAATATCAGTATCCACAGAAGGAGCAA  
 ATACTTCTCATCTACATCTACATCCACCCTGGGACAAGCCATCTTGTAAAATGTGCGGAGAAGGAGAA  
 AACTTTCTGTGTGAATGGAGGGGAGTGCTTCATGGTGAAGACCTTTCAAACCCCTCGAGATACTTGTGC  
 AAGTGCCCAAATGAGTTTACTGGTGATCGCTGCCAAAACCTACGTAATGGCCAGCTTCTACAAGGCGGAGG  
 AGCTGTACCAGAAGAGAGTGTGACCATAACCGGCATCTGCATCGCCCTCCTTGTGGTCCGCATCATGTG  
 TGTGGTGGCTACTGCAAAACCAAGAAACAGCGGAAAAAGCTGCATGACCGTCTTCGGCAGAGCCTTCGG  
 TCTGAACGAAACAATATGATGAACATTGCCAATGGGCCTCACCATCTAACCCACCCCGAGAATGTCC  
 AGCTGGTGAATCAATACGATCTAAAAACGTCATCTCCAGTGAGCATATTGTTGAGAGAGAAGCAGAGAC  
 ATCCTTTTCCACCAGTCACTATACTTCCACAGCCCATCACTCCACTACTGTCACCCAGACTCCTAGCCAC  
 AGCTGGAGCAACGACACACTGAAAGCATCCTTCCGAAAGCCACTCTGTAATCGTGATGCATCCGTAG  
 AAAACAGTAGGCACAGCAGCCCAACTGGGGGCCAAGAGGACGTCTTAATGGCACAGGAGGCCCTCGTGA  
 ATGTAACAGCTTCTCAGGCATGCCAGAGAAACCCTGATTCTACCGAGACTCTCCTCATAGTGAAGG  
 TATGTGTAGCCATGACCACCCGGCTCGTATGTACCTGTAGATTTCCACACGCCAAGCTCCCCAAAT  
 CGCCCCCTTCGGAATGTCTCCACCCGTGTCCAGCATGACGGTGTCCATGCCTTCCATGGCGGTGAGCCC  
 CTTTATGGAAGAAGAGAGACCTCTACTTCTCGTGACACCACCAAGGCTGCGGGAGAAGAAGTTTGACCAT  
 CACCCTCAGCAGTTCAGCTCCTTCCACCACAACCCGCGCATGACAGTAACAGCCTCCCTGCTAGCCCT  
 TGAGGATAGTGGAGGATGAGGAGTATGAAACGACCCAAGAGTACGAGCCAGCCCAAGAGCCTGTTAAGAA  
 ACTCGCCAATAGCCGGCGGGCCAAAAGAACCAAGCCCAATGGCCACATTGCTAACAGATTGGAAGTGGAC  
 AGCAACACAAGCTCCAGAGCAGTAACCTCAGAGAGTGAACAGAAAGATGAAAGAGTAGGTGAAGATACGC  
 CTTTCCTGGGCATACAGAACCCCTGGCAGCCAGTCTTGAGGCAACACCTGCCTTCCGCCTGGCTGACAG  
 CAGGACTAACCCAGCAGGCCGCTTCTCGACACAGGAAGAAATCCAGGCCAGGCTGTCTAGTGAATTGCT  
 AACCAAGACCCTATTGCTGTA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC220213 representing NM\_013957  
 Red=Cloning site Green=Tags(s)


MSEKKEGRGKGGKGGKKERGSQSPALPPRLKEMKSEQESAAGSKLVLRCESTSEYSSSLRFKW  
 FKNGNELNRKNKPQNIKIQQKPGKSELRINKASLADSGEYMKVISKLGNSANITIVESNEIITGMP  
 ASTEGAYVSSSPIRISVSTEGANTSSTSTSTGTSHLVKCAEKEKTCVNGGECFMVKDLSNPSRYLC  
 KCPNEFTGDRSQNYVMASFYKAEELYQKRVLITIGICIALLVGIMCVVAYCKTKKQRKCLHDLRQSLR  
 SERNNMMNIANGPHHPNPPENVQLVNQYVSKNVISSEHIVEREAETSFSSTSHYTSTAHHSTTVTQTSPH  
 SWSNGHTESILSESHSVIVMSSVENSRRHSSPTGGPRGRLNGTGGPRECNSFLRHARETPDSYRDSPHSER  
 YYSAMTTPARMSPVDFHTPSSPKSPPSEMSPPVSSMTVSMPSMAVSPFMEEERPLLVTAPPRLREKFKDH  
 HPQQFSSFHHNP AHDNSLPA SPLRIVEDEEYETTQEYEPAQEPVKLANSSRAKRTKPNGHIANRLEVD  
 SNTSSQSSNSESETEDERVGEDTPFLGIQNPLAASLEATPAFRLADSRTPAGRFSTQEEIQARLSSVIA  
 NQDPIAV

TRTRP L E Q K L I S E E D L A A N D I L D Y K D D D D K V

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk8070\\_g09.zip](https://cdn.origene.com/chromatograms/mk8070_g09.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_013957

**ORF Size:** 1911 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_013957.5](#)

**RefSeq Size:** 2490 bp

**RefSeq ORF:** 1914 bp

**Locus ID:** 3084

**UniProt ID:** [Q02297](#)

**Cytogenetics:** 8p12

**Domains:** Neuregulin, ig, IGc2, IG, EGF, EGF


**Protein Families:** Druggable Genome, Secreted Protein, Transcription Factors, Transmembrane

**Protein Pathways:** ErbB signaling pathway

**MW:** 70 kDa

**Gene Summary:** The protein encoded by this gene is a membrane glycoprotein that mediates cell-cell signaling and plays a critical role in the growth and development of multiple organ systems. An extraordinary variety of different isoforms are produced from this gene through alternative promoter usage and splicing. These isoforms are expressed in a tissue-specific manner and differ significantly in their structure, and are classified as types I, II, III, IV, V and VI. Dysregulation of this gene has been linked to diseases such as cancer, schizophrenia, and bipolar disorder (BPD). [provided by RefSeq, Apr 2016]

## Product images:



Circular map for RC220213