

Product datasheet for **RC220249L3V**

TTC9 (NM_015351) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TTC9 (NM_015351) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TTC9
Synonyms:	TTC9A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_015351
ORF Size:	666 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220249).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015351.1
RefSeq Size:	5217 bp
RefSeq ORF:	669 bp
Locus ID:	23508
UniProt ID:	Q92623
Cytogenetics:	14q24.2
MW:	24.2 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein that contains three tetratricopeptide repeats. The gene has been shown to be hormonally regulated in breast cancer cells and may play a role in cancer cell invasion and metastasis. [provided by RefSeq, Mar 2009]