

Product datasheet for **RC220262A1V**

Human OR7D4 (NM_001005191) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human OR7D4 (NM_001005191) AAV Particle
Tag:	Myc-DDK
Symbol:	OR7D4
Synonyms:	hg105; OR7D4P; OR19-7; OR19-B; OR19B
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC220262 representing NM_001005191 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAAGCAGAAAACCTTACAGAATTATCAAAATTTCTCCTCCTGGGACTCTCAGATGATCCTGAACTGC
AGCCCGTCCTCTTTGGGCTGTTCTGTCCATGTACCTGGTCACGGTGCTGGGGAACCTGCTCATCATTCT
GGCCGTCAGCTCTGACTCCACCTCCACACCCCATGTACTTCTCCTCTCCAACCTGTCCTTTGTTGAC
ATCTGTTTCATCTCCACCACAGTCCCAAGATGCTAGTGAGCATCCAGGCACGGAGCAAGACATCTCCT
ACATGGGGTGCCTCACTCAGGTGATTTTTTAATGATGTTTCTGGAATGGATACTTTCTACTGGCCGT
GATGGCCTATGACCGGTTTGTGGCCATCTGCCACCCACTGCACTACACGGTCATCATGAACCCCTGCCTC
TGTGGCCTCCTGGTTCTGGCATCTTGGTTCATCATTTTCTGGTCTCCCTGGTTCATATTCTACTGATGA
AGAGGTTGACCTTCTCCACAGGCACTGAGATTCGCATTTCTTCTGTGAACCGGCTCAGGTCTCAAGGT
GGCCTGCTCTAACACCCCTCTCAATAACATTGTCTTGATGTGGCCACGGCACTGCTGGGTGTGTTTCT
GTAGCTGGGATCCTCTTCTCTACTCTCAGATTGTCTCCTCCTTAATGGGAATGCTCTCCACCAAGGGCA
AGTACAAAGCCTTTTCCACCTGTGGATCTCACCTCTGTGTGGTCTCCTTGTCTATGGAACAGGACTTGG
GGTCTATCTGAGTTCTGCTGTGACCCATTCTTCCAGAGCAGCTCCACCGCCTCAGTGATGTACGCCATG
GTCACCCCATGCTGAACCCCTTCATCTACAGCCTGAGGAACAAGGATGTGAAGGGGGCCCTGGAAGAC
TCCTCAGCAGGGCCGACTCTTGTTCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC220262 representing NM_001005191 Red=Cloning site Green=Tags(s) MEAENLTLSKFLLLGLSDDPELQPVLFGLFLSMYLVTVLGNLLIILAVSSDSDLHTPMYFFLSNLSFVD ICFISTTVPKMLVSIQARSKDISYMGCLTQVYFLMMFAGMDTFLAVMAYDRFVAICHPLHYTVIMNPCL CGLLVLASWFIIFWFSLVHILLMKRLTFSTGTEIPHFCEPAQVLKVACSNLNNIVLYVATALLGVFP VAGILFSYSQIVSSLMGMSSTKGKYKAFSTCGSHLCVVSIFYGTGLGVYLSSAVTHSSQSSSTASVMYAM VTPMLNPFIFYSLRNKDYKGALERLLSRADSCP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001005191
ORF Size:	936 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001005191.1, NP_001005191.1</u>
RefSeq Size:	939 bp
RefSeq ORF:	939 bp
Locus ID:	125958
UniProt ID:	<u>Q8NG98</u>
Cytogenetics:	19p13.2
MW:	34.3 kDa