

Product datasheet for RC220272L2

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

ADCY3 (NM_004036) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: ADCY3 (NM_004036) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: ADCY3

Synonyms: AC-III; AC3; BMIQ19

Mammalian Cell None

Selection:

Vector: pLenti-C-mGFP (PS100071)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC220272).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_004036

ORF Size: 3432 bp



ADCY3 (NM_004036) Human Tagged Lenti ORF Clone - RC220272L2

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 004036.2</u>

RefSeq Size:4410 bpRefSeq ORF:3435 bp

UniProt ID: O60266

Cytogenetics: 2p23.3

Domains: CYCc

Locus ID:

Protein Families: Druggable Genome, Transmembrane

109

Protein Pathways: Calcium signaling pathway, Chemokine signaling pathway, Dilated cardiomyopathy, Gap

junction, GnRH signaling pathway, Melanogenesis, Olfactory transduction, Oocyte meiosis, Progesterone-mediated oocyte maturation, Purine metabolism, Vascular smooth muscle

contraction, Vibrio cholerae infection


MW: 129 kDa

Gene Summary: This gene encodes adenylyl cyclase 3 which is a membrane-associated enzyme and catalyzes

the formation of the secondary messenger cyclic adenosine monophosphate (cAMP). This protein appears to be widely expressed in various human tissues and may be involved in a number of physiological and pathophysiological metabolic processes. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2016]



Product images:



Circular map for RC220272L2