

## Product datasheet for **RC220298L4V**

### **BS69 (ZMYND11) (NM\_006624) Human Tagged ORF Clone Lentiviral Particle**

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	BS69 (ZMYND11) (NM_006624) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BS69
Synonyms:	BRAM1; BS69; MRD30
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006624
ORF Size:	1806 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220298).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006624.5</a>
RefSeq Size:	4336 bp
RefSeq ORF:	1809 bp
Locus ID:	10771
UniProt ID:	<a href="#">Q15326</a>
Cytogenetics:	10p15.3
Domains:	PWWP, BROMO, PHD, zf-MYND
Protein Families:	Druggable Genome, Transcription Factors


[View online »](#)

**MW:** 71 kDa

**Gene Summary:** The protein encoded by this gene was first identified by its ability to bind the adenovirus E1A protein. The protein localizes to the nucleus. It functions as a transcriptional repressor, and expression of E1A inhibits this repression. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]