

Product datasheet for **RC220303L4V**

PSD3 (NM_015310) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PSD3 (NM_015310) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PSD3
Synonyms:	EFA6D; EFA6R; HCA67
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015310
ORF Size:	3141 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220303).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015310.3
RefSeq Size:	11687 bp
RefSeq ORF:	3144 bp
Locus ID:	23362
UniProt ID:	Q9NYI0
Cytogenetics:	8p22
Domains:	Sec7, PH
Protein Pathways:	Endocytosis



[View online »](#)

MW: 115.8 kDa

Gene Summary: Guanine nucleotide exchange factor for ARF6.[UniProtKB/Swiss-Prot Function]