

## Product datasheet for **RC220307L3V**

### SLC39A12 (NM\_152725) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC39A12 (NM_152725) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC39A12
Synonyms:	bA570F3.1; LZT-Hs8; ZIP-12
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_152725
ORF Size:	1962 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220307).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_152725.1</a>
RefSeq Size:	2636 bp
RefSeq ORF:	1965 bp
Locus ID:	221074
UniProt ID:	<a href="#">Q504Y0</a>
Cytogenetics:	10p12.33
Domains:	Zip
Protein Families:	Transmembrane



[View online »](#)

**MW:** 72.7 kDa

**Gene Summary:** Zinc is an essential cofactor for hundreds of enzymes. It is involved in protein, nucleic acid, carbohydrate, and lipid metabolism, as well as in the control of gene transcription, growth, development, and differentiation. SLC39A12 belongs to a subfamily of proteins that show structural characteristics of zinc transporters (Taylor and Nicholson, 2003 [PubMed 12659941]).[supplied by OMIM, Aug 2008]