

## Product datasheet for **RC220354L3V**

### PSD3 (NM\_206909) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PSD3 (NM_206909) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PSD3
Synonyms:	EFA6D; EFA6R; HCA67
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_206909
ORF Size:	1539 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220354).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_206909.1</a> , <a href="#">NP_996792.1</a>
RefSeq Size:	6822 bp
RefSeq ORF:	1542 bp
Locus ID:	23362
UniProt ID:	<a href="#">Q9NYI0</a>
Cytogenetics:	8p22
Protein Pathways:	Endocytosis
MW:	57.9 kDa



[View online »](#)

**Gene Summary:** Guanine nucleotide exchange factor for ARF6.[UniProtKB/Swiss-Prot Function]