

Product datasheet for **RC220356L1V**

C12orf48 (PARPBP) (NM_017915) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | C12orf48 (PARPBP) (NM_017915) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | C12orf48 |
| Synonyms: | AROM; C12orf48; PARI |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_017915 |
| ORF Size: | 1737 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC220356). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_017915.1 |
| RefSeq Size: | 2175 bp |
| RefSeq ORF: | 1740 bp |
| Locus ID: | 55010 |
| UniProt ID: | Q9NWS1 |
| Cytogenetics: | 12q23.2 |
| Protein Families: | Druggable Genome |
| MW: | 65.05 kDa |



[View online »](#)

Gene Summary:

Required to suppress inappropriate homologous recombination, thereby playing a central role DNA repair and in the maintenance of genomic stability. Antagonizes homologous recombination by interfering with the formation of the RAD51-DNA homologous recombination structure. Binds single-strand DNA and poly(A) homopolymers. Positively regulate the poly(ADP-ribosyl)ation activity of PARP1; however such function may be indirect. [UniProtKB/Swiss-Prot Function]