

Product datasheet for **RC220359A1V**

Human CD79A (NM_001783) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CD79A (NM_001783) AAV Particle
Tag: Myc-DDK
Symbol: CD79A
Synonyms: IGA; MB-1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC220359 representing NM_001783
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCCTGGGGTCCAGGAGTCTCCAAGCTCTGCCTGCCACCATCTTCTCTCTTCTGCTGTCTGCTG
TCTACCTGGGCCCTGGGTGCCAGGCCCTGTGGATGCACAAGGTCCCAGCATATTGATGGTGAGCCTGGG
GGAAGACGCCCCTTCCAATGCCCGACAATAGCAGCAACAACGCCAACGTACCTGGTGGCGCGTCTCT
CATGGCAACTACACGTGGCCCCCTGAGTTCTTGGGCCCGGGCGAGGACCCCAATGGTACGCTGATCATCC
AGAATGTGAACAAGAGCCATGGGGGCATATACGTGTGCCGGTCCAGGAGGGCAACGAGTCATACCAGCA
GTCTGTGGGCACCTACCTCCGCGTGGCCAGCCGCCCGCCAGGCCCTTCTGGACATGGGGGAGGGCACC
AAGAACCGAATCATCACAGCCGAGGGGATCATCTCTGTTCTGCGCGGTGGTGCCTGGGACGCTGCTGC
TGTTTCAGGAAACGATGGCAGAACGAGAAGCTCGGGTTGGATGCCGGGGATGAATATGAAGATGAAAACCT
TTATGAAGGCCTGAACCTGGACGACTGCTCCATGTATGAGGACATCTCCCGGGGCTCCAGGGCACCTAC
CAGGATGTGGGCAGCCTCAACATAGGAGATGTCCAGCTGGAGAAGCCG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC220359 representing NM_001783
Red=Cloning site Green=Tags(s)

MPGGPGVLQALPATIFLLFLLSAVYLGPGQALWMHKVPASLMVSLGEDAHFQCPHNSSNNANVTWVRVL
HGNYTWPPEFLGPGEDPNGTLIIQNVNKGSHGGIYVCRVQEGNESYQQSCGTYLVRQPPPRPFLDMGEGT
KNRIITAEGIIILLFCAVVPGLTLLFRKRWQNEKLGDLADGEYEDENLYEGLNLDDCSMYEDI SRGLQGT
QDVGSLNIGDVQLEKP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_001783
ORF Size:	678 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001783.3
RefSeq Size:	1258 bp
RefSeq ORF:	681 bp
Locus ID:	973
UniProt ID:	P11912
Cytogenetics:	19q13.2
MW:	25.04 kDa