

Product datasheet for RC220365

OR11G2 (NM_001005503) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OR11G2 (NM_001005503) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	OR11G2
Synonyms:	OR14-34
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC220365 representing NM_001005503 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCATTTTCTTTCCCAAATGATTTAAATATAAATCTGATTCCCATCTATGTTTGACCAGTCAATCAG
TAATTGCTGGTGTCTTTACAATTCACAGGCACATGAAAATCTTCAACAGCCCCAGCAACTCCAGCACCT
CACTGGCTTCATCCTCCTGGGCTCCCTTGCCCCAGGGAGGGGCAGATCCTCCTTTGTGCTCTCACT
GTTGTTTACCTCCTGACCCTCATGGCAATGGTCCATCATCTGTGCTGTGCACTGGGATCAGAGACTCC
ACGCCCCATGTACATCCTGCTCGCCAATTCTCCTTCTGGAGATATGTTATGCACCTCCACAGTCCC
CAGCATGCTGGCCAACTCCTCCTGACACCAAGATCATCTCGTTCTCTGGCTGCTTCTCCAGTTCTAC
TTTTTCTTCTCCTTGGGCTCTACAGAATGCTTTTTCTGGCAGTTATGGCATTGATCGATACCTTGCCA
TCTGTGCGGCTCTACGCTATCCAACATTATGACCAGACGTCTCTGTACCAATCTGTGGTCAATTGCTG
GGTACTTGGTTTCATCTGGTTCTTGATTCTATCGTCAACATCTCCCAAATGTCTTCTGTGGATCAGG
ATTATTGACCACTTCTATGTGACCCAGCTCCTCTTCTAACTCTCACTTGCAAAAAGGCCCTGTGATAG
AGCTTGTCTTTCTGTCTTAAGTCCCTGCTGCTTTATGCTCTTTCTTTCATTGTGGGGTCCATGCT
TCTGGTGTGAGAGCTGTGTTGAGGGTCCCTTCAGCAGCTGGGAGAAGAAAGGCTTTCTCCACCTGTGG
TCTCACCTGGCTGTGGTTCACTGTTCTACGGCTCAGTACTGGTCATGTATGGGAGCCACCATCTAAGA
ATGAAGCTGAAAAGCAGAAGACTGTGACTCTGTTTTATTCTGTTGTTACCCCACTGTTAACCCCTGTGAT
ATATAGTCTTAGGAACAAAGATATGAGAAAAGCTCTGAAGAAATTTGGGGAACA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC220365 representing NM_001005503
Red=Cloning site Green=Tags(s)

MHFLSQNDLNINLIPHLCLHRHSVIAGAFTHRHMKIFNSPNSNSTFTGFILLGFPCPREGQILLFVLF
 VVYLLTLMGNGSIICAVHWDQRLHAPMYILLANFSFLEICYVTSTVPSMLANFLSDTKIISFSGCFLQFY
 FFFSLGSTECFFLAVMAFDRYLAICRPLRYPTIMTRRLCTNLVNCWVLGFIWFLIPIVNI SQMSFCGSR
 IIDHFLCDPAPLLTLTCKKGPVIELVFSVLSPLPVFMLFLFIVGSYALVVRVAVLRVPSAAGRRAKAFSTCG
 SHLAVVSLFYGSVLMYGSPPSKNEAGKQKTVTLFYSVVTPLLNPVIYSLRNKDMRKALKKFWGT


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8004_b02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



ACCN: NM_001005503

ORF Size: 1035 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001005503.1](#), [NP_001005503.1](#)

RefSeq Size: 1038 bp

RefSeq ORF: 936 bp

Locus ID: 390439

UniProt ID: [Q8NGC1](#)

Cytogenetics: 14q11.2


Protein Families: Transmembrane

Protein Pathways: Olfactory transduction

MW: 38.7 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC220365