

Product datasheet for **RC220369**

GNLY (NM_012483) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: GNLY (NM_012483) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: GNLY
Synonyms: D2S69E; LAG-2; LAG2; NKG5; TLA519
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC220369 representing NM_012483
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAAGGTCTGGTCTTCTCTCGTCTGAGCCCTGAGTACTACGACCTGGCAAGAGCCCACCTGCGTGATG
AGGAGAAATCCTGCCGTGCCTGGCCAGGAGGGCCCCAGGGTGACCTGTTGACCAAACACAGGAGCT
GGGCCGTGACTACAGGACCTGTCTGACGATAGTCCAAAACCTGAAGAAGATGGTGGATAAGCCCACCCAG
AGAAGTGTTCGAATGCTGCGACCCGGTGTGTAGGACGGGAGGTACGATGGCGGACGTCTGCAGAA
ATTCATGAGGAGGTATCAGTCTAGAGTACCCAGGGCCTCGTGGCCGGAGAACTGCCAGCAGATCTG
TGAGGACCTCAGGTTGTATACCTTCTACAGGTCCCCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC220369 representing NM_012483
Red=Cloning site **Green**=Tags(s)

MEGLVFSRLSPEYYDLARAHLRDEEKSCPCLAQEGPQGDLTKTQELGRDYRTCLTIVQKLKMKMVDKPTQ
RSVSNAATRVCRTRGRSRWRDVCNFMRRYQSRVTQGLVAGETAQQICEDLRLCIPSTGPL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_012483

ORF Size: 390 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_012483.4](#)

RefSeq Size: 995 bp

RefSeq ORF: 393 bp

Locus ID: 10578

UniProt ID: [P22749](#)


Cytogenetics: 2p11.2

Protein Families: Secreted Protein

MW: 14.7 kDa

Gene Summary: The product of this gene is a member of the saposin-like protein (SAPLIP) family and is located in the cytotoxic granules of T cells, which are released upon antigen stimulation. This protein is present in cytotoxic granules of cytotoxic T lymphocytes and natural killer cells, and it has antimicrobial activity against *M. tuberculosis* and other organisms. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC220369