

## Product datasheet for RC220393L2

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### Cytokeratin 18 (KRT18) (NM\_000224) Human Tagged Lenti ORF Clone

**Product data:** 

**Product Type:** Expression Plasmids

**Product Name:** Cytokeratin 18 (KRT18) (NM\_000224) Human Tagged Lenti ORF Clone

Tag: mGFP

Symbol: Cytokeratin 18

Synonyms: CK-18; CYK18; K18

Mammalian Cell None

Selection:

Vector:pLenti-C-mGFP (PS100071)E. coli Selection:Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(RC220393).

Sequence:

**Restriction Sites:** Sgfl-Mlul

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF.

ACCN: NM\_000224

ORF Size: 1290 bp



#### Cytokeratin 18 (KRT18) (NM\_000224) Human Tagged Lenti ORF Clone - RC220393L2

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 000224.2</u>

 RefSeq Size:
 1485 bp

 RefSeq ORF:
 1293 bp

 Locus ID:
 3875

 UniProt ID:
 P05783

Cytogenetics: 12q13.13

**Domains:** filament

**Protein Pathways:** Pathogenic Escherichia coli infection

**MW:** 47.9 kDa

**Gene Summary:** KRT18 encodes the type I intermediate filament chain keratin 18. Keratin 18, together with its

filament partner keratin 8, are perhaps the most commonly found members of the


intermediate filament gene family. They are expressed in single layer epithelial tissues of the

body. Mutations in this gene have been linked to cryptogenic cirrhosis. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul

2008]



# **Product images:**



Circular map for RC220393L2