

Product datasheet for RC220393L4V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Cytokeratin 18 (KRT18) (NM_000224) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: Cytokeratin 18 (KRT18) (NM 000224) Human Tagged ORF Clone Lentiviral Particle

Symbol: Cytokeratin 18

Synonyms: CK-18; CYK18; K18

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

ACCN: NM_000224 **ORF Size:** 1290 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC220393).

Sequence:
OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 000224.2

RefSeq Size: 1485 bp
RefSeq ORF: 1293 bp
Locus ID: 3875
UniProt ID: P05783
Cytogenetics: 12q13.13
Domains: filament

Protein Pathways: Pathogenic Escherichia coli infection





MW: 47.9 kDa

Gene Summary: KRT18 encodes the type I intermediate filament chain keratin 18. Keratin 18, together with its

filament partner keratin 8, are perhaps the most commonly found members of the intermediate filament gene family. They are expressed in single layer epithelial tissues of the body. Mutations in this gene have been linked to cryptogenic cirrhosis. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul

2008]