

## Product datasheet for **RC220434L2V**

### Calcipressin 1 (RCAN1) (NM\_004414) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Calcipressin 1 (RCAN1) (NM_004414) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Calcipressin 1
Synonyms:	ADAPT78; CSP1; DSC1; DSCR1; MCIP1; RCN1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_004414
ORF Size:	756 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220434).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004414.5</a>
RefSeq Size:	2457 bp
RefSeq ORF:	759 bp
Locus ID:	1827
UniProt ID:	<a href="#">P53805</a>
Cytogenetics:	21q22.12
Domains:	Calcipressin
Protein Families:	Transcription Factors



[View online »](#)

**MW:** 27.9 kDa

**Gene Summary:** The protein encoded by this gene interacts with calcineurin A and inhibits calcineurin-dependent signaling pathways, possibly affecting central nervous system development. This gene is located in the minimal candidate region for the Down syndrome phenotype, and is overexpressed in the brain of Down syndrome fetuses. Chronic overexpression of this gene may lead to neurofibrillary tangles such as those associated with Alzheimer disease. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2013]