

## Product datasheet for **RC220443L3V**

### **RNASE4 (NM\_002937) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	RNASE4 (NM_002937) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RNASE4
Synonyms:	RAB1; RNS4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002937
ORF Size:	441 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220443).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_002937.3</a>
RefSeq Size:	1808 bp
RefSeq ORF:	444 bp
Locus ID:	6038
UniProt ID:	<a href="#">P34096</a>
Cytogenetics:	14q11.2
Domains:	RNAse_Pc
Protein Families:	Secreted Protein, Transmembrane



[View online »](#)

**MW:** 16.8 kDa

**Gene Summary:** The protein encoded by this gene belongs to the pancreatic ribonuclease family. It plays an important role in mRNA cleavage and has marked specificity towards the 3' side of uridine nucleotides. Alternative splicing results in four transcript variants encoding the same protein. This gene and the gene that encodes angiogenin share promoters and 5' exons. Each gene splices to a unique downstream exon that contains its complete coding region. [provided by RefSeq, Aug 2013]