

Product datasheet for **RC220455L1V**

RPL38 (NM_000999) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RPL38 (NM_000999) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RPL38
Synonyms:	L38
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_000999
ORF Size:	210 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220455).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000999.2
RefSeq Size:	368 bp
RefSeq ORF:	213 bp
Locus ID:	6169
UniProt ID:	P63173
Cytogenetics:	17q25.1
Domains:	Ribosomal_L38e
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Ribosome

MW: 8 kDa

Gene Summary: Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L38E family of ribosomal proteins. It is located in the cytoplasm. Alternative splice variants have been identified, both encoding the same protein. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome, including one located in the promoter region of the type 1 angiotensin II receptor gene. [provided by RefSeq, Jul 2008]