

## Product datasheet for **RC220487L4V**

### MAPKAP Kinase 2 (MAPKAPK2) (NM\_004759) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	MAPKAP Kinase 2 (MAPKAPK2) (NM_004759) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MAPKAP Kinase 2
Synonyms:	MAPKAP-K2; MK-2; MK2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004759
ORF Size:	1110 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220487).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004759.3</a>
RefSeq Size:	3608 bp
RefSeq ORF:	1113 bp
Locus ID:	9261
UniProt ID:	<a href="#">P49137</a>
Cytogenetics:	1q32.1
Domains:	pkinese, TyrKc, S_TKc
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

**Protein Pathways:** MAPK signaling pathway, Neurotrophin signaling pathway, VEGF signaling pathway

**MW:** 42 kDa

**Gene Summary:** This gene encodes a member of the Ser/Thr protein kinase family. This kinase is regulated through direct phosphorylation by p38 MAP kinase. In conjunction with p38 MAP kinase, this kinase is known to be involved in many cellular processes including stress and inflammatory responses, nuclear export, gene expression regulation and cell proliferation. Heat shock protein HSP27 was shown to be one of the substrates of this kinase in vivo. Two transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]