

Product datasheet for **RC220517A1V**

Human NM23A (NME1) (NM_198175) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human NM23A (NME1) (NM_198175) AAV Particle
Tag:	Myc-DDK
Symbol:	NM23A
Synonyms:	AWD; GAAD; NB; NBS; NDKA; NDPK-A; NDPKA; NM23; NM23-H1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC220517 representing NM_198175 Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC ATGGTGCTACTGTCTACTTTAGGGATCGTCTTTCAAGGCGAGGGGCCTCCTATCTCAAGCTGTGATACAGGAACCATGGCCAACTGTGAGCGTACCTTCATTGCGATCAAACCAGATGGGGTCCAGCGGGGTCTTGTTGGGAGAGATTATCAAGCGTTTTGAGCAGAAAGGATTCGCCCTTGTTGGTCTGAAATTCATGCAAGCTTCCGAAGATCTTCTCAAGGAACACTACGTTGACCTGAAGGACCGTCCATTCTTTGCCGGCCTGGTGAATATACATGCCTCAGGCGCGGTAGTTGCCATGGTCTGGGAGGGGCTGAATGTGGTGAAGACGGGCCGAGTCATGCTCGGGAGACCAACCCTGCAGACTCCAAGCCTGGGACCATCCGTGGAGACTTCTGCATACAAGTTGGCAGGAACATTATACATGGCAGTGATTCTGTGGAGAGTGCAGAGAAGGAGATCGGCTTGTGGTTTCACCTGAGGAACGGTAGATTACACGAGCTGTGCTCAGAACTGGATCTATGAA ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>RC220517 representing NM_198175 Red=Cloning site Green=Tags(s) MVLLSTLGIVFQGEPPISSCDTGTMANCERTFIAIKPDGVQRGLVGEI IKRFEQKGFRLVGLKFMQASEDLLKEHYVDLDRPFFAGLVKYMHS GPVVMVWEGLNVVKTGRVMLGETNPADSKPGTIRGDFCIQVGRNIIHGSDSVESA EKEIGLWFHPEELVDYTSQAQNWIE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2


[View online »](#)

ACCN:	NM_198175
ORF Size:	531 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_198175.1</u> , <u>NP_937818.1</u>
RefSeq Size:	1031 bp
RefSeq ORF:	534 bp
Locus ID:	4830
UniProt ID:	<u>P15531</u>
Cytogenetics:	17q21.33
MW:	19.5 kDa