

Product datasheet for **RC220528L3V**

SIDT2 (NM_001040455) Human Tagged ORF Clone Lentiviral Particle

Product data:


Product Type:	Lentiviral Particles
Product Name:	SIDT2 (NM_001040455) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SIDT2
Synonyms:	CGI-40
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001040455
ORF Size:	2496 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220528).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001040455.1
RefSeq Size:	3919 bp
RefSeq ORF:	2499 bp
Locus ID:	51092
UniProt ID:	Q8NBJ9
Cytogenetics:	11q23.3
Protein Families:	Transmembrane
MW:	94.3 kDa



[View online »](#)

Gene Summary:

Mediates the translocation of RNA and DNA across the lysosomal membrane during RNA and DNA autophagy (RDA), a process in which RNA or DNA is directly imported into lysosomes in an ATP-dependent manner, and degraded (PubMed:27046251, PubMed:27846365). Involved in the uptake of single-stranded oligonucleotides by living cells, a process called gymnosis (PubMed:28277980). Involved in the uptake of single-stranded oligonucleotides by living cells, a process called gymnosis. In vitro, mediates the uptake of linear DNA more efficiently than that of circular DNA, but exhibits similar uptake efficacy toward RNA and DNA. Binds long double-stranded RNA (dsRNA) (500 - 700 base pairs), but not dsRNA shorter than 100 bp (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:

[RC220528L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC220528L3V particle to overexpress human SIDT2-Myc-DDK fusion protein.