

Product datasheet for **RC220544A1V**

Human Cyclin C (CCNC) (NM_005190) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Cyclin C (CCNC) (NM_005190) AAV Particle
Tag:	Myc-DDK
Symbol:	Cyclin C
Synonyms:	CycC; hSRB11; SRB11
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC220544 representing NM_005190 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCAGGGAACCTTTGGCAGAGCTCCCACTATTTGCAATGGATTTGGATAAACAAGATCTGTTGAAGG
AGCGCCAAAAGGATTTAAAGTTTCTCTCAGAGGAAGAATATTGGAAGTTACAAATATTTTTACAAATGT
TATCCAAGCATTAGGTGAACATCTTAAATTAAGACAACAAGTTATTGCCACTGCTACGGTATATTTCAAG
AGATTCTATGCCAGGTATTCTCTGAAAAGTATAGATCCTGTATTAATGGCTCCTACATGTGTGTTTTGG
CATCCAAAGTAGAGGAATTTGGAGTAGTTTCAAATACAAGATTGATTGCTGCTGCTACTTCTGTATTA
AACTAGATTTTTCATATGCCTTTCCAAAGGAATTTCTTATAGGATGAATCATATATTAGAATGTGAATTC
TATCTGTTAGAACTAATGGATTGTTGCTTGATAGTGTATCATCCTTATAGACCTTTGCTCCAGTATGTGC
AGGACATGGGCCAAGAAGACATGTTGCTTCCCCTTGCATGGAGGATAGTGAATGATACCTACAGAACGGA
TCTTTGCCTACTGTATCCTCCTTTCATGATAGCTTTAGCTTGCCTACATGTAGCCTGTGTTGTACAGCAG
AAAGATGCCAGGCAATGGTTTGCTGAGCTTTCTGTGGATATGAAAAGATTTGGAAATAATCAGGGTTA
TTTTAAACTATATGAGCAGTGAAGAATTTTCATGAGAGAAAAGAGATGGCAACCATTCTAGTAAGAT
GCCAAACCAAAACACCTCCAAACAGTGAAGGAGAGCAGGGTCCAAATGGAAGTCAGAACTCTAGCTAC
AGCCAATCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC220544 representing NM_005190 Red=Cloning site Green=Tags(s) MAGNFWQSSHYLQWILDKQDLLKERQKDLKFLSEEEYWKLQIFFTNVIQALGEHLKLRQQVIATATVYFK RFYARYSLKSIDPVLMAPTCVFLASKVEEFGVVSNTLIAAATSVLKTRFSYAFPKFPPYRMNHILECEF YLLELMDCCLIYVHPYRPLLQYVQDMGQEDMLLPLAWRIVNDTYRTDLCLLYPPFMIALACLVACVVQQ KDARQWFAELSVDMEKILEIIRVILKLYEQWKNFDERKEMATILSKMPKPKPPNSEGEQGPNGSQNSSY SQS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_005190
ORF Size:	849 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_005190.3</u>
RefSeq Size:	2321 bp
RefSeq ORF:	852 bp
Locus ID:	892
UniProt ID:	<u>P24863</u>
Cytogenetics:	6q16.2
MW:	33.1 kDa