

Product datasheet for **RC220560L4V**

GPR31 (NM_005299) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GPR31 (NM_005299) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GPR31
Synonyms:	12-HETER; HETER; HETER1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005299
ORF Size:	957 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220560).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005299.1
RefSeq Size:	960 bp
RefSeq ORF:	960 bp
Locus ID:	2853
UniProt ID:	O00270
Cytogenetics:	6q27
Protein Families:	Druggable Genome, GPCR, Transmembrane
MW:	34.9 kDa



[View online »](#)

Gene Summary:

High-affinity receptor for 12-(S)-hydroxy-5,8,10,14-eicosatetraenoic acid (12-S-HETE). 12-(S)-HETE is an arachidonic acid metabolite secreted by platelets and tumor cells, and known to induce endothelial cells retraction allowing invasive cell access to the subendothelial matrix, which is a critical step for extravasation or metastasis. Ligand-binding lead to activation of ERK1/2 (MAPK3/MAPK1), MEK, and NF-kappa-B.[UniProtKB/Swiss-Prot Function]