

Product datasheet for **RC220572L2V**

AK3L1 (AK4) (NM_013410) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AK3L1 (AK4) (NM_013410) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AK3L1
Synonyms:	AK3; AK3L1; AK3L2; AK 4
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_013410
ORF Size:	669 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220572).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013410.2
RefSeq Size:	2199 bp
RefSeq ORF:	672 bp
Locus ID:	205
UniProt ID:	P27144
Cytogenetics:	1p31.3
Domains:	ADK, ADK_lid
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Metabolic pathways, Purine metabolism

MW: 25.1 kDa

Gene Summary: This gene encodes a member of the adenylate kinase family of enzymes. The encoded protein is localized to the mitochondrial matrix. Adenylate kinases regulate the adenine and guanine nucleotide compositions within a cell by catalyzing the reversible transfer of phosphate group among these nucleotides. Five isozymes of adenylate kinase have been identified in vertebrates. Expression of these isozymes is tissue-specific and developmentally regulated. A pseudogene for this gene has been located on chromosome 17. Three transcript variants encoding the same protein have been identified for this gene. Sequence alignment suggests that the gene defined by NM_013410, NM_203464, and NM_001005353 is located on chromosome 1. [provided by RefSeq, Jul 2008]