

## Product datasheet for **RC220581L3V**

### FTCD (NM\_206965) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	FTCD (NM_206965) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FTCD
Synonyms:	LCHC1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_206965
ORF Size:	1623 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220581).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_206965.1</a>
RefSeq Size:	1884 bp
RefSeq ORF:	1626 bp
Locus ID:	10841
UniProt ID:	<a href="#">O95954</a>
Cytogenetics:	21q22.3
Protein Pathways:	Histidine metabolism, Metabolic pathways, One carbon pool by folate
MW:	58.7 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a bifunctional enzyme that channels 1-carbon units from formiminoglutamate, a metabolite of the histidine degradation pathway, to the folate pool. Mutations in this gene are associated with glutamate formiminotransferase deficiency. Alternatively spliced transcript variants have been found for this gene.[provided by RefSeq, Dec 2009]