

Product datasheet for **RC220613L3V**

PHC2 (NM_198040) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PHC2 (NM_198040) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PHC2
Synonyms:	EDR2; HPH2; PH2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_198040
ORF Size:	2574 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220613).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_198040.2 , NP_932157.1
RefSeq Size:	3873 bp
RefSeq ORF:	2577 bp
Locus ID:	1912
UniProt ID:	Q8IXK0
Cytogenetics:	1p35.1
MW:	90.7 kDa



[View online »](#)

Gene Summary:

In *Drosophila melanogaster*, the 'Polycomb' group (PcG) of genes are part of a cellular memory system that is responsible for the stable inheritance of gene activity. PcG proteins form a large multimeric, chromatin-associated protein complex. The protein encoded by this gene has homology to the *Drosophila* PcG protein 'polyhomeotic' (Ph) and is known to heterodimerize with EDR1 and colocalize with BMI1 in interphase nuclei of human cells. The specific function in human cells has not yet been determined. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]