

Product datasheet for **RC220631L3V**

ZNF274 (NM_016325) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF274 (NM_016325) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF274
Synonyms:	HFB101; ZF2; ZKSCAN19; ZSCAN51
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016325
ORF Size:	1863 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220631).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016325.2 , NP_057409.1
RefSeq Size:	2743 bp
RefSeq ORF:	1866 bp
Locus ID:	10782
UniProt ID:	Q96GC6
Cytogenetics:	19q13.43
Protein Families:	Transcription Factors
Protein Pathways:	Neurotrophin signaling pathway



[View online »](#)

MW: 70.3 kDa

Gene Summary: This gene encodes a zinc finger protein containing five C2H2-type zinc finger domains, one or two Kruppel-associated box A (KRAB A) domains, and a leucine-rich domain. The encoded protein has been suggested to be a transcriptional repressor. It localizes predominantly to the nucleolus. Alternatively spliced transcript variants encoding different isoforms exist. These variants utilize alternative polyadenylation signals. [provided by RefSeq, Jul 2008]