

## Product datasheet for **RC220669L1V**

### **TNFAIP8 (NM\_001077654) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	TNFAIP8 (NM_001077654) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TNFAIP8
Synonyms:	GG2-1; MDC-3.13; NDED; SCC-S2; SCCS2
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001077654
ORF Size:	564 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220669).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001077654.1</a>
RefSeq Size:	1979 bp
RefSeq ORF:	567 bp
Locus ID:	25816
UniProt ID:	<a href="#">O95379</a>
Cytogenetics:	5q23.1
Protein Families:	Druggable Genome
MW:	21.7 kDa



[View online »](#)

**Gene Summary:**

Acts as a negative mediator of apoptosis and may play a role in tumor progression. Suppresses the TNF-mediated apoptosis by inhibiting caspase-8 activity but not the processing of procaspase-8, subsequently resulting in inhibition of BID cleavage and caspase-3 activation.[UniProtKB/Swiss-Prot Function]