

## Product datasheet for **RC220673L2V**

### **BTN3A3 (NM\_006994) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	BTN3A3 (NM_006994) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BTN3A3
Synonyms:	BTF3; BTN3.3
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_006994
ORF Size:	1752 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220673).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006994.3</a>
RefSeq Size:	2956 bp
RefSeq ORF:	1755 bp
Locus ID:	10384
UniProt ID:	<a href="#">O00478</a>
Cytogenetics:	6p22.2
Domains:	IGv, IG, SPRY, PRY
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

**MW:** 64.8 kDa

**Gene Summary:** The butyrophilin (BTN) genes are a group of major histocompatibility complex (MHC)-associated genes that encode type I membrane proteins with 2 extracellular immunoglobulin (Ig) domains and an intracellular B30.2 (PRYSPRY) domain. Three subfamilies of human BTN genes are located in the MHC class I region: the single-copy BTN1A1 gene (MIM 601610) and the BTN2 (e.g., BTN2A1; MIM 613590) and BTN3 (e.g., BNT3A3) genes, which have undergone tandem duplication, resulting in 3 copies of each (summary by Smith et al., 2010 [PubMed 20208008]).[supplied by OMIM, Nov 2010]