

## Product datasheet for **RC220708L4V**

### Zinc finger protein 287 (ZNF287) (NM\_020653) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Zinc finger protein 287 (ZNF287) (NM_020653) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Zinc finger protein 287
Synonyms:	ZKSCAN13; ZSCAN45
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020653
ORF Size:	2262 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220708).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020653.1</a> , <a href="#">NP_065704.1</a>
RefSeq Size:	4278 bp
RefSeq ORF:	2286 bp
Locus ID:	57336
UniProt ID:	<a href="#">Q9HBT7</a>
Cytogenetics:	17p11.2
Domains:	KRAB, LER, zf-C2H2
Protein Families:	Transcription Factors



[View online »](#)

**MW:** 87.6 kDa

**Gene Summary:** This gene encodes a member of the krueppel family of zinc finger proteins, suggesting a role as a transcription factor. Its specific function has not been determined. This gene is located near the Smith-Magenis syndrome region on chromosome 17. [provided by RefSeq, Jul 2008]