

Product datasheet for **RC220754A1V**

Human IL17E (IL25) (NM_172314) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human IL17E (IL25) (NM_172314) AAV Particle

Tag: Myc-DDK

Symbol: IL17E

Synonyms: IL17E

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC220754 representing NM_172314
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTACCAGGTGGTTGCATTCTTGGCAATGGTCATGGGAACCCACACCTACAGCCACTGGCCCAGCTGCT
GCCCCAGCAAAGGGCAGGACACCTCTGAGGAGCTGCTGAGGTGGAGCACTGTGCCTGTGCCTCCCTAGA
GCCTGTAGGCCCAACCGCCACCCAGAGTCCTGTAGGGCCAGTGAAGATGGACCCCTCAACAGCAGGGCC
ATCTCCCCCTGGAGATATGAGTTGGACAGAGACTTGAACCGGCTCCCCAGGACCTGTACCACGCCGTT
GCCTGTGCCCGCACTGCGTCAGCCTACAGACAGGCTCCACATGGACCCCGGGGCAACTCGGAGCTGCT
CTACCACAACCAGACTGTCTTCTACCGGCGGCCATGCCATGGCGAGAAGGGCACCCACAAGGGCTACTGC
CTGGAGCGCAGGCTGTACCGTGTTCCTTAGCTTGTGTGTGTGTGCGGCCCGTGTGATGGGC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC220754 representing NM_172314
Red=Cloning site Green=Tags(s)

MYQVVAFLAMVMGTHYSHWPSCCPKQDTSEELLRWSTVPVPPLEPARPNRHPESCRASEDGPLNSRA
ISPWRYELDRDLNRLPQDLYHARCLCPHCVSLQTGSHMDPRGNSELLYHNQTVFYRRPCHGEKGTHTKGYC
LERRLYRVSLACVCVRPRVMG

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2



[View online »](#)

ACCN:	NM_172314
ORF Size:	483 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_172314.1, NP_758525.1</u>
RefSeq Size:	1187 bp
RefSeq ORF:	486 bp
Locus ID:	64806
UniProt ID:	<u>Q9H293</u>
Cytogenetics:	14q11.2
MW:	16.7 kDa