

## Product datasheet for **RC220756L3V**

### TRAR4 (TAAR6) (NM\_175067) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TRAR4 (TAAR6) (NM_175067) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TAAR6
Synonyms:	TA4; taR-4; taR-6; TAR4; TAR6; TRAR4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_175067
ORF Size:	1035 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220756).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_175067.1</a> , <a href="#">NP_778237.1</a>
RefSeq Size:	1038 bp
RefSeq ORF:	1038 bp
Locus ID:	319100
UniProt ID:	<a href="#">Q96RI8</a>
Cytogenetics:	6q23.2
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction



[View online »](#)

**MW:** 38.3 kDa

**Gene Summary:** This gene encodes a seven-transmembrane G-protein-coupled receptor that likely functions as a receptor for endogenous trace amines. Mutations in this gene may be associated with schizophrenia.[provided by RefSeq, Feb 2010]