

Product datasheet for RC220780L3

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

NF-kB p65 (RELA) (NM_021975) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: NF-kB p65 (RELA) (NM_021975) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: NF-kB p65

Synonyms: CMCU; NFKB3; p65

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Sgfl-Mlul


E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC220780).

Sequence:

Restriction Sites:
Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_021975

ORF Size: 1653 bp

NF-kB p65 (RELA) (NM_021975) Human Tagged Lenti ORF Clone - RC220780L3

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 021975.2</u>

 RefSeq Size:
 1760 bp

 RefSeq ORF:
 1656 bp

 Locus ID:
 5970

 UniProt ID:
 004206

Cytogenetics: 11q13.1

Domains: RHD. IPT

Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: Acute myeloid leukemia, Adipocytokine signaling pathway, Apoptosis, B cell receptor signaling

pathway, Chemokine signaling pathway, Chronic myeloid leukemia, Cytosolic DNA-sensing pathway, Epithelial cell signaling in Helicobacter pylori infection, MAPK signaling pathway, Neurotrophin signaling pathway, NOD-like receptor signaling pathway, Pancreatic cancer, Pathways in cancer, Prostate cancer, RIG-I-like receptor signaling pathway, Small cell lung

cancer, T cell receptor signaling pathway, Toll-like receptor signaling pathway

MW: 60 kDa


Gene Summary: NF-kappa-B is a ubiquitous transcription factor involved in several biological processes. It is

held in the cytoplasm in an inactive state by specific inhibitors. Upon degradation of the inhibitor, NF-kappa-B moves to the nucleus and activates transcription of specific genes. NF-kappa-B is composed of NFKB1 or NFKB2 bound to either REL, RELA, or RELB. The most abundant form of NF-kappa-B is NFKB1 complexed with the product of this gene, RELA. Four transcript variants encoding different isoforms have been found for this gene. [provided by

RefSeq, Sep 2011]



Product images:



Circular map for RC220780L3