

Product datasheet for **RC220821L4V**

NAP1L1 (NM_004537) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NAP1L1 (NM_004537) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NAP1L1
Synonyms:	NAP1; NAP1L; NRP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004537
ORF Size:	1173 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220821).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004537.3 , NP_004528.1
RefSeq Size:	2908 bp
RefSeq ORF:	1176 bp
Locus ID:	4673
UniProt ID:	P55209
Cytogenetics:	12q21.2
Domains:	NAP
MW:	45.2 kDa


[View online »](#)

Gene Summary:

This gene encodes a member of the nucleosome assembly protein (NAP) family. This protein participates in DNA replication and may play a role in modulating chromatin formation and contribute to the regulation of cell proliferation. Alternative splicing results in multiple transcript variants encoding different isoforms; however, not all have been fully described. [provided by RefSeq, Apr 2015]