

Product datasheet for **RC220823**

ZYG11B (NM_024646) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZYG11B (NM_024646) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZYG11B
Synonyms:	ZYG11
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC220823 representing NM_024646
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGCCCGAGGACCAGGCCGGCGCAGCCATGGAGGAGCGTCTCCCTATTCCTTACTTGATATCTGCTTGA
 ATTTCTTGACTACTCACCTTGAGAAGTTCTGTTACGCCAGACAAGATGGAACATTGTGTCTGCAGGAACC
 TGGAGTATTCCACAGGAGGTGGCTGATCGACTGCTTCGACCATGGCTTTTCATGGTCTATTGAATGAT
 GGAAGTGTGGGTATTTTAGGGGCAACCAGATGCGCTTAAAGCGAGCCTGCATTGCGAAAGCAAAGATCT
 CTGCTGTTGCTTCCGGAAAGCTTTCTGCCACCACAAGTTAGTGAACCTTGATGCCACAGGTGTGAATGC
 TGATATCACGATTACAGACATTATCAGTGGGCTTGGCAGTAACAAATGGATCCAGCAGAATCTCCAGTGC
 CTGGTGTGAATTCATTAACCTCTCTCCCTCGAGGATCCTACGAGCGCTGCTTCAGCCGGCTTTCTGGCC
 TTCGAGCTTTAAGCATCACGAATGTTCTCTTTTACAATGAAGACCTGGCTGAAGTTGCCTCATTGCCAAG
 ATTAGAGAGCTTGGATATTTCTAACACCTCAATCACAGACATCACTGCTCTACTGGCCTGCAAAGACCGA
 CTCAAGTCTCTAACCATGCACCATTGAAATGTTTAAAAATGACAACATCCAGATACTGGATGTAGTTC
 GGAACCTCAAACATCTGAATCATCTTGATATCTCAGATGATAAACAGTTTACATCAGACATAGCTTTTCG
 CTTACTAGAACAAAAAGACATCCTACCTAACCTTGTCTCTGGATGTTTCTGGGAGAAAGCACGTGACA
 GATAAAGCCGTTGAAGCCTTTATACAACAACGTCCAAGCATGCAATTTGTAGGTTTGTGGCTACTGATG
 CTGGTTACTCTGAATTCCTCACAGGCGAAGGACATTTGAAGGTGTCTGGGAAGCCAATGAAACTCAGAT
 TGCAGAAGCACTGAAGCGTTACAGTGAACGGGCATTCTTTGTTGCGGAAGCTCTATTTTCATCTTTTGT
 CTGACTCATGTGATGGAAAAACAAAGCCAGAAATTTTAAAGCTTGTGGTTACTGGGATGAGAAACCACC
 CTATGAATTTGCCAGTGAACCTGGCTGCAAGCGCCTGTGATTTAACTTAACCAAGCAGGATCTTGCTGC
 AGGGATGCCTGTCCGACTCCTGGCTGATGTGACCCATTTGCTGCTCAAAGCCATGGAACATTTTCCCAAT
 CACCAGCAGTTACAGAAGAATTGCCTCCTTTCACTTTGCAGTGACCGGATCCTTCAAGATGTTCCATTTA
 ACAGGTTTGAAGCAGCCAAGCTTGTGATGAGTGGCTTTGCAACCATGAGGATCAAAACATGCAAAGGAT
 GGCAGTTGCTATCATTTCTATCCTGGCTGCCAAGCTTTCTACAGAACAACTGCACAGCTTGGTACTGAG
 CTCTTCATTGTCAGGCAACTTCTTCAAATAGTGAAGCAGAAAACCAATCAAAATTCAGTGGACACTACAT
 TGAAATTTACTTTGAGTGCATTTGGAACCTCACAGATGAATCTCCAACCACTTGTAGACACTTTATTGA
 AAACCAAGGGTTAGAATCTTCATGAGGGTTCTAGAGTCTTTCCCACTGAGTCATCCATTCAGCAGAAA
 GTTCTAGGACTTTTGAACAATATAGCTGAAGTACAAGAATTACATTCTGAATTAATGTGGAAAGATTTTA
 TAGACCACATCAGTAGTCTCCTACACAGTGTGGAAGTGAAGTCAGTTACTTTGCAGCTGGAATTATTGC
 CCATTTAATATCCAGAGGTGAACAAGCTTGGACATTGAGTCGTAGCCAGAGGAATTCTCTGCTGGATGAT
 TTGCATTACAGCTATTTTGAAATGGCCAACTCCAGAGTGTGAGATGGTAGCATACAGGTCCTTTAATCCAT
 TTTTCCCATTAATTGGCTGTTTCAACACCAGGAGTTTCAGCTATGGGCAGTTTGGGCCATGCAACATGT
 CTGCAGCAAGAATCCTTCAAGGTATTGCAGCATGCTGATTGAAGAAGGAGGATTGCAGCATTTATACAAC
 ATCAAAGATCATGAACATACTGATCCCCATGTCCAACAGATTGCTGTGGCCATTCTGGATAGCTTAGAAA
 AACACATTGTGCGCCATGGGAGGCCACCTCCCTGTAAAAAACAGCCCCAAGCCAGACTAAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC220823 representing NM_024646
 Red=Cloning site Green=Tags(s)


MPEDQAGAAMEEASPYSLLDICLNFLTTHLEKFC SARQDGTLC LQEPGVFPQEVADRLRLTMAFHGLLND
 GTVGIFRGNQMRLKRACIRKAKISAVAFKAFCHHKLVELDATGVNADITITDIIISGLGSNKWIQNLQC
 LVLNSLTLSLEDPYERCF SRLSGLRAL SITNVLFYNEDLAEVASLPRLSLDISNTSITDITALLACKDR
 LKSLTMHHLKCLKMTTQTILDVVRELKHLNHLDISDDKQFTSDIALRLLEQKDILPNLVSLDVSGRKHVT
 DKAVEAFIQQRPSMQFVGLLATDAGYSEFLTGEHLKVSGEANETQIAEALKRYSERAFFVREALFHLFS
 LTHVMEKTKPEILKL VVTGMRNHPMNL PVQLAASACVFNLTKQDLAAGMPVRLLADVTHLLLKAMEHFPN
 HQQLQKNCLLSLCSDRILQDVPFNRFEAAKLVQWLCNHEDQNMQRMAVAIISILA AKLSTEQTAQLGTE
 LFIVRQLLQIVKQKTNQNSVDTLKFTLSALWNLDESPTTCRHF IENQGLELFMRVLESFPTESSIQQK
 VLGLLNIAEVQELHSELMWKDFIDHISLLHSVEVEVSFYAAGIIAHLISRGEQAWTLSRSQRNSLLDD
 LHSAILKWPTPECEMVAYRSFNPFPLLCFTTPGVQLWAVWAMQHVCSKNPSRYCSMLIEEGGLQHLYN
 IKDHEHTDPHVQQIAVAILDSEKHI VRHGRPPPCKKQPQARLN

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8041_d11.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_024646


ORF Size: 2232 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_024646.3</u>
RefSeq Size:	6404 bp
RefSeq ORF:	2235 bp
Locus ID:	79699
UniProt ID:	<u>Q9C0D3</u>
Cytogenetics:	1p32.3
MW:	83.7 kDa
Gene Summary:	Probably acts as target recruitment subunit in the E3 ubiquitin ligase complex ZYG11B-CUL2-Elongin BC.[UniProtKB/Swiss-Prot Function]


Product images:



Circular map for RC220823



Western blot validation of overexpression lysate (Cat# [LY411165]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC220823 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified ZYG11B protein (Cat# [TP320823]). The protein was produced from HEK293T cells transfected with ZYG11B cDNA clone (Cat# RC220823) using MegaTran 2.0 (Cat# [TT210002]).