

Product datasheet for **RC220823L4V**

ZYG11B (NM_024646) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZYG11B (NM_024646) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZYG11B
Synonyms:	ZYG11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024646
ORF Size:	2232 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC220823).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024646.1
RefSeq Size:	6404 bp
RefSeq ORF:	2235 bp
Locus ID:	79699
UniProt ID:	Q9C0D3
Cytogenetics:	1p32.3
MW:	83.7 kDa


[View online »](#)

Gene Summary:

Probably acts as target recruitment subunit in the E3 ubiquitin ligase complex ZYG11B-CUL2-Elongin BC.[UniProtKB/Swiss-Prot Function]